

S3



The S3 Scanner



What It Is

The patented Pharmanex BioPhotonic Scanner is a cutting edge testing tool that noninvasively measures carotenoid levels in living tissue.

The S3 Scanner is fast, small, truly portable, and wirelessly controlled with an Apple device.



What It Does

It provides an immediate indication of a person's overall carotenoid levels. It helps you track your Skin Carotenoid Score (SCS) over time as you make important diet and lifestyle changes.

Scan anywhere, anytime.



How It Works

Using optical signals obtained through Resonance Raman Spectroscopy, you can get a reading of the carotenoid levels in your skin.

Receive your SCS score by placing your hand in front of the scanner's blue light. See your score in 30 seconds.



Your Score

LOW CAROTENOID ABSORPTION INTO TISSUE

Red 10,000-19,000
Orange 20,000-29,000

AVERAGE CAROTENOID ABSORPTION INTO TISSUE

Yellow 30,000-39,000
Green 40,000-49,000

ABOVE AVERAGE CAROTENOID ABSORPTION INTO TISSUE

Blue 50,000-59,000
Platinum 60,000-89,000+



Increase Your Score

YOUR SCORE INCREASES AS YOU

- Consume more fruits and vegetables
- Manage a healthy and low-stress lifestyle



Learn more about the S3 Scanner at nuskin.com - Opportunity - Scanner

